


POSTERIOR SURGICAL LIQUIDS



**“Unmatched purity
of Vitreco
surgical liquids.”**

Vitreco a company of:



BVI



Vitreco

Perfluorocarbons


Vitrex developed two perfluorocarbons for the vitreoretinal surgery. Deca-Vit and Octa-Vit for draining sub retinal fluids and flattening retina in vitreoretinal surgery. Additional indications include: lifting diving nucleus or luxated IOL's, bleeding control etc.



Packaged box/5, Sterile

- Naturally-pure basic materials, additionally highly purified.
- Maximum interfacial tension to avoid mixture with other liquids.
- Syringe packaged, ready-for-use.
- 5 ml and 7 ml syringe packaging.
- Produced according to EN ISO 13485-2012.
- CE certified.

Deca-Vit




Specifications

- Specific Gravity: 1.93 G/cm³
- Refractive Index: 1.31
- Boiling Point: 140 Celsius

Surgical Benefits

- Higher gravity
- High(er) interfacial tension to silicone oil
- High(er) refractive index for visible differentiation with BSS and silicone oil

Octa-Vit



Specifications

- Specific Gravity: 1.77 G/cm³
- Refractive Index: 1.27
- Boiling Point: 101-105 Celsius

Surgical Benefits

- High vapor pressure allows evaporation of residual droplets
- Easier complete removal at the end of the surgery

Silicone oils


Vitreoq developed for the vitreoretinal surgery two silicone oils Sil-Vit 1.000 and Sil-Vit 5.000 for use as an endo-tamponade to re-attach detached retina to the underlying choroid.

- Naturally-pure basic materials, additionally highly purified.
- Maximum interfacial tension to avoid mixture with other liquids.
- Syringe packaged, ready-for-use.
- Produced according to EN ISO 13485-2012.
- CE certified.



Packaged box/5, Sterile

Sil-Vit 1000

syringe	<p>Art.No.: SL00.S11</p>  <p>Sil-Vit 1000 Cst Syringe, 10 ml</p>
---------	---


Specifications

- Viscosity: 1,000 – 1,500 mPas (Centistokes)
- Refractive Index: 1.40
- Specific Gravity: 0.97 G/cm³
- Surface Tension: 21 mn/m versus air and 40 mn/m versus water

Surgical Benefits

- Low viscosity, easy and quick to inject and to remove

Sil-Vit 5000

syringe	<p>Art.No.: SL00.S13</p>  <p>Sil-Vit 5000 Cst Syringe, 10 ml</p>
---------	--

Specifications

- Viscosity: 5,000 – 5,900 mPas (Centistokes)
- Refractive Index: 1.40
- Specific Gravity: 0.97 G/cm³
- Surface Tension: 21 mn/m versus air and 40 mn/m versus water

Surgical Benefits

- 5,000 centistokes viscosity for high viscosity endo-tamponade



ISO Certification

Vitreu's quality system complies with the international quality and regulatory requirements and is certified by the Notified Body DQS Medizinprodukte GmbH, according to international standard ISO 13485 Medical devices – Quality management systems – Requirements for regulatory purposes.

CE Certification

The Vitrea medical devices are certified according to the medical device directive 93/42/EEC.

European Countries

To order Vitrea products in Europe contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Ltd

85B Park Drive, Milton Park
Abingdon, Oxfordshire
OX14 4RY, United Kingdom
E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Belgium

E: bestellingen@bvimedical.com
T: 32.52.280518 *Option 3*
F: 32.52.890243

France, Luxembourg

E: commandes@bvimedical.com
T: 33.4.85880033 *Option 3*
F: 33.4.88714057

Germany, Switzerland, Austria

E: bestellungen@bvimedical.com
T: 49.622.164.79854 *Option 3*
F: 49.622.139.29883

Italy

E: ordini@bvimedical.com
T: 39.02.89877520 *Option 3*
F: 39.02.44386868

Netherlands

E: bestellingen@bvimedical.com
T: 31.20808.3202 *Option 3*
F: 31.20203.1293

Spain

E: ordenes@bvimedical.com
T: 34.91.1230335 *Option 3*
F: 34.91.1849957

United Kingdom, Eire

E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

All Other Regions

E: ordersrow@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Vitreu Headquarters

Seggelant-Noord 2
3237 MG Vierpolders
The Netherlands
E: info@vitrea.com
T: 31.181.745090
F: 31.181.478583



Vitreu a company of:



vitrea.com bvimedical.com

Authorized Vitrea distributor:

