

KUHN 45° POWER FORCEPS

Patent Pending



“I have not come across any task in the eye that I am not able to do with the Kuhn 45° Power Forceps.”

Professor Ferenc Kuhn

The Kuhn 45° Power Forceps is designed for peeling proliferative intravitreal, epi and sub-retinal membranes, lens capsule, iris, small IOFB, fibrinous membranes and more. The forceps is also suitable for peeling standard ILM and ERM membranes and attached posterior hyaloids.

Its unique tip design enables exceptional visualization while grasping the target tissue, with a 45° crocodile design inner surface that enables selective and strong grasping of the target tissue, regardless of its rigidity, size, orientation, or surface texture.

Available in 23G and 25G in both reusable (R06) and disposable (D06) versions. The reusable forceps will be delivered including the unique cleaning and sterilization holder. The disposable versions are packaged box/5, Sterile.

23 Gauge

Art.No.: FC23.R06
Art.No.: FC23.D06



23G Forceps Kuhn 45°
Power Gripping


25 Gauge

Art.No.: FC25.R06
Art.No.: FC25.D06



25G Forceps Kuhn 45°
Power Gripping

Vitrex a company of:



BVI



Vitrex



ISO Certification

Vitre^q's quality system complies with the international quality and regulatory requirements and is certified by the Notified Body DQS Medizinprodukte GmbH, according to international standard ISO 13485 Medical devices – Quality management systems – Requirements for regulatory purposes.

CE Certification

The Vitre^q medical devices are certified according to the medical device directive 93/42/EEC.

FDA

Vitre^q B.V. is an FDA registered manufacturer of medical devices.

European Countries

To order Vitre^q products in Europe contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Ltd

85B Park Drive, Milton Park
Abingdon, Oxfordshire
OX14 4RY, United Kingdom
E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Belgium

E: bestellingen@bvimedical.com
T: 32.52.280518 *Option 3*
F: 32.52.890243

France, Luxembourg

E: commandes@bvimedical.com
T: 33.4.85880033 *Option 3*
F: 33.4.88714057

Germany, Switzerland, Austria

E: bestellingen@bvimedical.com
T: 49.622.16479854 *Option 3*
F: 49.622.13929883

Italy

E: ordini@bvimedical.com
T: 39.02.89877520 *Option 3*
F: 39.02.44386868

Netherlands

E: bestellingen@bvimedical.com
T: 31.20808.3202 *Option 3*
F: 31.20203.1293

Spain

E: ordenes@bvimedical.com
T: 34.91.1230335 *Option 3*
F: 34.91.1849957

United Kingdom, Eire

E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

North and South America

Ordering Vitre^q products is easy. Contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Inc.

500 Totten Pond Rd – 10 CityPoint
Waltham, MA 02451 USA

United States, Canada

E: customersupport@bvimedical.com
T: 1.866.906.8080 *Toll-free*
781.906.8080 *Direct*
F: 1.866.906.4304

Latin America, all other

E: customersupport_lam@bvimedical.com
T: 781.906.8080 *Ext. 7*
F: 781.906.4304

Japan

E: bvi.cs@bvimedical.com
T: 03.5651.2860
F: 03.3662.2242

All Other Regions

E: ordersrow@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Vitre^q Headquarters

Seggelant-Noord 2
3237 MG Vierpolders
The Netherlands
E: info@vitreq.com
T: 31.181.745090
F: 31.181.478583



Vitre^q a company of:



vitreq.com bvimedical.com

Authorized Vitre^q distributor:

