




Atraumatic Retinal Surface Brushing and Manipulation

Retractable Tissue Manipulator

Patent Pending



- Rigid Silicone tip structure without diamond particles – eliminates any risk of losing particles inside the eye
- Retractable silicone tip ensures smooth passage through entry systems

Vitreco
Vitreoretinal Surgical Products

Retractable Tissue Manipulator

INDICATIONS FOR USE



The tissue manipulator can be used in cases of retinal disorders to manipulate the retina and retinal folds.

FEATURES

- Retractable silicone tip with special ribbed structure
- The green colored silicone tip allows for a clear differentiation to the red retinal surface

Tips

Available in 23G and 25G. Disposable and ready to use.

Reference	Article Name	Packaging	
TM23.D01	23G Retractable Tissue Manipulator	Packaged box/5, Sterile	
TM25.D01	25G Retractable Tissue Manipulator	Packaged box/5, Sterile	

“The more I work with it, the more I have the impression that it is an indispensable instrument.”*

— Professor Hassan Mortada



*Read the full interview on our website

FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT