



29G SPOTLIGHT DIRECTIONAL CHANDELIER


Patent Pending



Directional function

**“A bright vision
on illumination”**

Vitreco a company of:



BVI




Vitreco

29G Spotlight Directional Chandelier

A new generation chandelier with major improvements in functionality and ease of use.

Patent Pending

The Vitreo 29G Spotlight Directional Chandelier featuring high-end specifications with consistent performance for bimanual ophthalmic endo-illumination procedures.



Directional function




“The Spotlight Directional Chandelier offers an easy and quick insertion in combination with an effective directional control of the light beam.”

Developed in collaboration with Professor Claus Eckardt

Benefits

- Improved illumination / visualization due to directional illumination system.
- Directional Chandelier facilitates bimanual surgery.
- Adapters provide compatibility to most endo-illumination systems.



29G entry system



29G minimally invasive incision

Features

- Directional Chandelier system designed to enhance illumination control.
- Wide view fiber tip for global endo-illumination.
- Easy insertion and fixation of the fiber with 29G valved entry system.
- 29G minimal invasive incision.
- Integrated scleral marker




29G Spotlight Directional Chandelier

Art.No.: LF29.D02
Packaged box/5, Sterile



Light Source Adaptors



The Vitreo 29G Spotlight Directional Chandelier is compatible with most endo-illumination systems in combination with Light Source Adaptors:

- Constellation® Light Source Adaptors
- Stellaris® Light Source Adaptors
- Geuder® Light Source Adaptors
- Ruck® Light Source Adaptors
- DORC® Light Source Adaptors

<p>Art.No.: LF00.R01</p>  <p>Constellation® Light Source Adaptor</p>	<p>Art.No.: LF00.R02</p>  <p>Stellaris® Light Source Adaptor</p>	<p>Art.No.: LF00.R03</p>  <p>Geuder® Light Source Adaptor</p>
<p>Art.No.: LF00.R04</p>  <p>Ruck® Light Source Adaptor</p>	<p>Art.No.: LF00.R05</p>  <p>DORC® Light Source Adaptor</p>	

*All product names, trademarks and registered trademarks are property of their respective owners.

Endo-Illumination Brief Statement for Vitreo Spotlight Directional Chandelier:

Indications for Use: Vitreo Spotlight Chandelier is intended for endo-illumination during ophthalmic surgery.

Precautions: Do not use product if package integrity has been breached or compromised; do not use product after the expiration date.

Caution (USA): Federal law restricts this device to sale by or on order of a licensed physician.

Attention: See Instructions for Use included in product packaging for complete listing of indications and warnings.





ISO Certification

Vitreo's quality system complies with the international quality and regulatory requirements and is certified by the Notified Body DQS Medizinprodukte GmbH, according to international standard ISO 13485 Medical devices – Quality management systems – Requirements for regulatory purposes.

CE Certification

The Vitreo medical devices are certified according to the medical device directive 93/42/EEC.

FDA

Vitreo is an FDA registered manufacturer of medical devices.

European Countries

To order Vitreo products in Europe contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Ltd

85B Park Drive, Milton Park
Abingdon, Oxfordshire
OX14 4RY, United Kingdom
E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Belgium

E: bestellingen@bvimedical.com
T: 32.52.280518 *Option 3*
F: 32.52.890243

France, Luxembourg

E: commandes@bvimedical.com
T: 33.4.85880033 *Option 3*
F: 33.4.88714057

Germany, Switzerland, Austria

E: bestellingen@bvimedical.com
T: 49.622.16479854 *Option 3*
F: 49.622.13929883

Italy

E: ordini@bvimedical.com
T: 39.02.89877520 *Option 3*
F: 39.02.44386868

Netherlands

E: bestellingen@bvimedical.com
T: 31.20808.3202 *Option 3*
F: 31.20203.1293

Spain

E: ordenes@bvimedical.com
T: 34.91.1230335 *Option 3*
F: 34.91.1849957

United Kingdom, Eire

E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

North and South America

Ordering Vitreo products is easy. Contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Inc.

500 Totten Pond Rd – 10 CityPoint
Waltham, MA 02451 USA

United States, Canada

E: customersupport@bvimedical.com
T: 1.866.906.8080 *Toll-free*
781.906.8080 *Direct*
F: 1.866.906.4304

Latin America, all other

E: customersupport_lam@bvimedical.com
T: 781.906.8080 *Ext. 7*
F: 781.906.4304

Japan

E: bvi.cs@bvimedical.com
T: 03.5651.2860
F: 03.3662.2242

All Other Regions

E: ordersrow@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Vitreo Headquarters

Seggelant-Noord 2
3237 MG Vierpolders
The Netherlands
E: info@vitreo.com
T: 31.181.745090
F: 31.181.478583



Vitreo a company of:



vitreo.com bvimedical.com

Authorized Vitreo distributor:

