

SHAH X-TRA GRIP FORCEPS



“The Shah X-tra Grip Forceps employ tips that are fine enough to engage the edge of a membrane and the serrated platforms hold the tissue firmly in place.”

Gaurav Shah, MD, The Retina Institute

Dr. Shah has no financial interest with the product mentioned.

Features

- Fine Tipped Forceps with Serrated Jaws platform.
- Serrations angled 45° for maximum grasping of tissue.*

Benefits

- Fine tips allow for ILM lifting.
- Serrated platform grasps and holds tissue, preventing slippage.
- Effective for all types of tissue manipulation, ILM, ERM, PVR and diabetic tissue.

* Patent pending, Ferenc Kuhn, MD

Available in 23G, 25G and 27G, in both reusable (R07) and disposable (D07) versions. Each reusable forceps is delivered with a unique sterilization/cleaning holder. The disposable versions are packaged box/5, Sterile.

23 Gauge

Art.No.: FC23.R07
Art.No.: FC23.D07



23G Forceps
Shah X-tra Grip

25 Gauge

Art.No.: FC25.R07
Art.No.: FC25.D07



25G Forceps
Shah X-tra Grip


27 Gauge

Art.No.: FC27.R07
Art.No.: FC27.D07



27G Forceps
Shah X-tra Grip

Vitrex a company of:



BVI



Vitrex



ISO Certification

Vitre^q's quality system complies with the international quality and regulatory requirements and is certified by the Notified Body DQS Medizinprodukte GmbH, according to international standard ISO 13485 Medical devices – Quality management systems – Requirements for regulatory purposes.

CE Certification

The Vitre^q medical devices are certified according to the medical device directive 93/42/EEC.

FDA

Vitre^q B.V. is an FDA registered manufacturer of medical devices.

European Countries

To order Vitre^q products in Europe contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Ltd

85B Park Drive, Milton Park
Abingdon, Oxfordshire
OX14 4RY, United Kingdom
E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Belgium

E: bestellingen@bvimedical.com
T: 32.52.280518 *Option 3*
F: 32.52.890243

France, Luxembourg

E: commandes@bvimedical.com
T: 33.4.85880033 *Option 3*
F: 33.4.88714057

Germany, Switzerland, Austria

E: bestellungen@bvimedical.com
T: 49.622.164.79854 *Option 3*
F: 49.622.139.29883

Italy

E: ordini@bvimedical.com
T: 39.02.89877520 *Option 3*
F: 39.02.44386868

Netherlands

E: bestellingen@bvimedical.com
T: 31.20808.3202 *Option 3*
F: 31.20203.1293

Spain

E: ordenes@bvimedical.com
T: 34.91.1230335 *Option 3*
F: 34.91.1849957

United Kingdom, Eire

E: ordersuk@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

North and South America

Ordering Vitre^q products is easy. Contact a Customer Service Representative via e-mail, telephone or fax.

Beaver-Visitec International, Inc.

500 Totten Pond Rd – 10 CityPoint
Waltham, MA 02451 USA

United States, Canada

E: customersupport@bvimedical.com
T: 1.866.906.8080 *Toll-free*
781.906.8080 *Direct*
F: 1.866.906.4304

Latin America, all other

E: customersupport_lam@bvimedical.com
T: 781.906.8080 *Ext. 7*
F: 781.906.4304

Japan

E: bvi.cs@bvimedical.com
T: 03.5651.2860
F: 03.3662.2242

All Other Regions

E: ordersrow@bvimedical.com
T: 44.1865.601.256 *Option 3*
F: 44.1865.595.761

Vitre^q Headquarters

Seggelant-Noord 2
3237 MG Vierpolders
The Netherlands
E: info@vitreq.com
T: 31.181.745090
F: 31.181.478583



Vitre^q a company of:



vitreq.com bvimedical.com

Authorized Vitre^q distributor:

